

This economical and disposable fecal collection device is designed for simple and hygienic collection of fecal material in clinic or by your client. The two-part plastic container is designed for use as both a fecal collection container and a flotation exam container. Various gastrointestinal parasite eggs can be detected through this simple test by the use of flotation solutions. Each dispensing box of 50 fecal diagnostic kits also contains one pad of 50 client instruction sheets. Easy-to-use laboratory instructions are printed on the box.

ITEM#	DESCRIPTION	SIZE
601026	Fecal Diagnostic Kits	Box of 50
602026	Fecal Diagnostic Kit Refill	Bag of 50



Rev. 02/18