

Rapid Differential Stain Kit

This dip quick, three-step rapid differential stain kit is for use in routine differential staining of blood smears, blood parasites, microbiological organisms and bone marrow specimens. The stains have extensive applicability in cytology for use in fine-needle aspiration analysis. The kit will process a specimen in 30 seconds and utilizes a fixative and two staining solutions to achieve the desired results. Staining intensity can be adjusted to suit your requirements by adjusting the times in one or both staining solutions.

ITEM#	DESCRIPTION	SIZE
608041	Rapid Differential Stain Kit	3 x 180 mL
603041	Stain Kit Refill (Fixative-Red-Blue)	3 x 500 mL
604041	Stain Fixative #1	500 mL
605041	Stain Red #2	500 mL
606041	Stain Blue #3	500 mL



Available exclusively through your veterinarian from:



(800) 824-3703
www.mvivet.com



(888) 722-2242
www.aahamarketlink.com