

Topical

VetOne® Surgical Adhesive contains purified n-Butyl Cyanoacrylate with a D&C-certified dye that allows for differentiation between the product and bodily fluid, but it dries clear for a clean appearance. It sets up quickly (5-15 seconds) to form a thin, flexible waterproof bandage that is sloughed naturally in 5–7 days. Disposable applicator tips (included) fit snugly into the tip on the bottle and allow greater control over the drop size and placement of the product.

For topical application only. Not intended for internal closures or for use on infected or deep puncture wounds.

ITEM#	DESCRIPTION	SIZE
510099	Surgical Adhesive	3 mL

Designed for wound management in situations such as:

- Lacerations and abrasions
- Cat declaws
- Cuts
- Tail docking

NEW!

- **Larger 3 mL Size**
- **Easy to Squeeze Bottle**

VetOne® Surgical Adhesive bottle is redesigned in a larger size to hold more glue and provide more dispensing comfort.

Vibrant Blue Color Identification

VetOne® Surgical Adhesive contains purified n-Butyl Cyanoacrylate with a D&C certified dye that allows for differentiation between the product and bodily fluid.



Fast Dry Time

VetOne® Surgical Adhesive dries clear for clean appearance. It sets up quickly (5-15 seconds) and dries flexible for maintaining closure on wound site.

Sterile & Accurate Applicators

The disposable applicators come in a sterile pack, which compared with the competition, can make a difference in an effort to reduce bacterial migration into the surgery site. The tip fits snugly into the tip on the bottle and allows greater control over the size and placement to the drop of product.

Rev. 10/18