

Endotracheal Tubes

VetOne endotracheal tubes are made of medical-grade clear PVC and fitted with 15 mm standard connector. The specially designed cuff, unusual bonding process, and soft tip help ensure smooth intubation. Packaged sterile; single-use. Size indicates inner diameter (ID).

- Murphy eye
- Pilot balloon with integral one-way valve
- High-volume, low-pressure cuff
- Full-length radiopaque line
- 1 cm graduation markings
- 2 cm indicator mark
- Can be used orally or nasally

ITEM#	DESCRIPTION	SIZE
601029	Endotracheal Tube	3.0 mm
602029	Endotracheal Tube	3.5 mm
603029	Endotracheal Tube	4.0 mm
604029	Endotracheal Tube	4.5 mm
605029	Endotracheal Tube	5.0 mm
606029	Endotracheal Tube	5.5 mm
607029	Endotracheal Tube	6.0 mm
608029	Endotracheal Tube	6.5 mm
609029	Endotracheal Tube	7.0 mm
610029	Endotracheal Tube	7.5 mm
611029	Endotracheal Tube	8.0 mm
612029	Endotracheal Tube	8.5 mm
613029	Endotracheal Tube	9.0 mm
614029	Endotracheal Tube	9.5 mm
615029	Endotracheal Tube	10.0 mm

Endotracheal Tube Size Comparison


All measurements are in millimeters.

VetOne

ITEM#	ID	OD	CUFF
601029	3.0	4.2	7.3
602029	3.5	4.9	8.9
603029	4.0	5.5	8.9
604029	4.5	6.3	11.0
605029	5.0	6.7	13.0
606029	5.5	7.3	16.0
607029	6.0	8.0	16.0
608029	6.5	8.7	20.0
609029	7.0	9.3	20.0
610029	7.5	10.0	24.0
611029	8.0	10.7	24.0
612029	8.5	11.3	26.0
613029	9.0	12.0	28.0
614029	9.5	12.7	28.0
615029	10.0	13.3	28.0

Sheridan

ITEM#	ID	OD	CUFF
000730	3.0	4.2	7.3
005944	3.5	4.9	8.9
005945	4.0	5.5	8.9
005946	4.5	6.2	11.0
005947	5.0	6.8	13.0
005948	5.5	7.5	16.0
005951	6.0	8.2	16.0
005952	6.5	8.9	19.0
005953	7.0	9.6	19.0
005954	7.5	10.2	24.0
005955	8.0	10.9	24.0
005956	8.5	11.5	27.0
005957	9.0	12.2	27.0
005958	9.5	12.9	26.0
007993	10.0	13.5	27.0



Available exclusively through your veterinarian from:



(800) 824-3703
www.mwivet.com



(888) 722-2242
www.aahamarketlink.com